

User Documentation

(HDAD: Heart Diagnostics and Analytics Dashboard)

System Requirements


You may run the system web interface application any OS including Windows XP, Windows Vista, Windows 7, MAC OS, Android, IOS, or Unix Firefox and Chrome

Access to the web application

Any user having internet connection and using one the listed browsers can access the web interface application using the following link: <https://hdad.herokuapp.com>

1. Homepage

- a. This is the first page displayed on the app.
- b. This Page contains some information about our app, a button to the sign-up page and some Navigation links at top of the page.



2. Sign Up Page (FR1-FR2)

- a. The access to the application is restricted to authorized users only so the user has to create an account to use the application.
- b. Enter the details in the input boxes.
- c. After submitting you will be redirected to welcome page and your username will be sent to email provided by you.
- d. Now you can go to Sign In page by clicking on navbar on the top of application.

[Home](#)[Sign In](#)[Sign Up](#)

Create Your Account

Name

Hospital

Email

Password

Confirm Password

Minimum 6 and Maximum 80 characters

[SIGN UP](#)

[Home](#)[Sign In](#)[Sign Up](#)

Welcome to HDAD

Your Username has been send to the Email provided by you.

3. Sign in Page (FR3, FR12)

- The access to the application is restricted to authorized users only, each user should know his username and password.
- You can Sign in by entering your username and password.
- If the username and password are correct you will be redirected to the dashboard.
- If the any or both of password and username is wrong you will see a message saying "Invalid Username or Password."

[Home](#)[Sign In](#)[Sign Up](#)

Sign In

Username

Password

☒ Remember me

SIGN IN

[Home](#)[Sign In](#)[Sign Up](#)

Sign In

Invalid username or password

Username

Password

☒ Remember me

SIGN IN

4. Dashboard (FR4, FR7)

- The user can perform many tasks on the dashboard.

[Home](#)[Dashboard](#)[Name](#)[Logout](#)

Quick Prediction

CHF DETECTION

HF CLASSIFICATION

HEART DISEASE

COMPLETE DIAGNOSTICS

Patient Record

S.no	Patient Id	Patient Name	Date Added	CHF	Heart Failure	Heart Disease
------	------------	--------------	------------	-----	---------------	---------------

Add Patient

Patient Name

ADD

BROWSE

SUBMIT

DOWNLOAD TEMPLATE

a. Adding Patient (FR5)

- User can add patients to the application by providing the name of the patient.
- The Patient Name will be visible in the table along with the results of the last diagnostics of the patient.

Home Dashboard Name Logout

Quick Prediction

CHF DETECTION HF CLASSIFICATION

HEART DISEASE COMPLETE DIAGNOSTICS

Patient Record

S.no	Patient Id	Patient Name	Date Added	CHF	Heart Failure	Heart Disease
------	------------	--------------	------------	-----	---------------	---------------

Home Dashboard Name Logout

Quick Prediction

CHF DETECTION HF CLASSIFICATION

HEART DISEASE COMPLETE DIAGNOSTICS




Patient Record

S.no	Patient Id	Patient Name	Date Added	CHF	Heart Failure	Heart Disease
1	PA0348	Patient Name	20/04/2021	-	-	-

b. Quick Predictions (FR4)

- User can run quick predictions by providing correct details on the respective forms.
- The result of the prediction will be visible to the user.


Note: - The result of these predictions will not be saved




c. Adding Patients in Bulk (Added Feature)

The User can add multiple patients at once by following these steps: -


- i. User can download the template for the file that has to be filled and uploaded by clicking on download template button.






- ii. A zip file will be downloaded that contain an example file and an instruction file to help fill the template properly.




- iii. After filling the csv file, the user can click on browse and select the csv and click on submit to start adding multiple patients at once.






d. Patient Diagnostics (FR6)

- i. User can click on the patient id button to see patient details and run diagnostics on that particular Patient.



5. Patient Details Page

- The User can perform many tasks on this page



a. Predictions (FR5, FR6)

- This is same as quick predictions the only difference is that here the user can save the data if he/she wants to.
- User can run quick predictions by providing correct details on the respective forms.
- The result of the prediction will be visible to the user.

Home

DashboardNameLogout

PA0445

Name : Vanessa Russell

Date Created : 20/04/2021

Time : 22:04:25

Diagnostics Date : 20/04/2021

NOTIFY PATIENT

Pulse (beats per minute)

80

Brain Natriuretic Peptide (pg/ml)

320

Cystatin (mg/l)

2

Potassium (mmol/L)

6.5

New York Heart Association (NYHA) Functional Classification

2

Killip Grade

1

SUBMIT

SAVE

ATIONHEART DISEASE

ESUPLOAD

CHF

Yes

Heart Failure

HFpEF

Heart Disease

No

Home

DashboardNameLogout

PA0445

Name : Vanessa Russell

Date Created : 20/04/2021

Time : 22:04:25

Diagnostics Date : 20/04/2021

NOTIFY PATIENT

Enter Details For CHF Detection

Patient Have Congestive Heart Failure

Pulse (beats per minute)

80

Brain Natriuretic Peptide (pg/ml)

320

Cystatin (mg/l)

2

Potassium (mmol/L)

6.5

New York Heart Association (NYHA) Functional Classification

2

ATIONHEART DISEASE

ESUPLOAD

CHF

Yes

Heart Failure

HFpEF

Heart Disease

No

Home

DashboardNameLogout

PA0445

Name : Vanessa Russell

Date Created : 20/04/2021

Time : 22:04:25

Diagnostics Date : 20/04/2021

NOTIFY PATIENT

Patient Have Congestive Heart Failure

Pulse (beats per minute)

80

Brain Natriuretic Peptide (pg/ml)

320

Cystatin (mg/l)

2

Potassium (mmol/L)

6.5

New York Heart Association (NYHA) Functional Classification

2

Killip Grade

1

SUBMIT

SAVE

ATIONHEART DISEASE

ESUPLOAD

CHF

Yes

Heart Failure

HFpEF

Heart Disease

No

Home
Dashboard
Name
Logout

PA0445

Name : Vanessa Russell
Date Created : 20/04/2021
Time : 22:04:25
Diagnostics Date : 20/04/2021

NOTIFY PATIENT

CHF DETECTION
HF CLASSIFICATION
HEART DISEASE

DIAGNOSTIC
STATISTICS
UPLOAD

S.no	Date	Time	CHF	Heart Failure	Heart Disease
1	20/04/2021	22:04:25	Yes	-	-
2	20/04/2021	22:04:25	Yes	HFpEF	No

b. Bulk Diagnostics (Added Feature)

The User can add multiple diagnostics at once by following these steps: -

- All the buttons for this function are available inside the popup that is shown when user click on upload.
- User can download the template for the file that has to be filled and uploaded by clicking on download template button.
- A zip file will be downloaded that contain an example file and an instruction file to help fill the template properly.
- After filling the csv file, the user can click on browse and select the csv and click on submit to start adding multiple diagnostics at once.

Home
Dashboard
Name
Logout

PA0445

Name : Vanessa Russell
Date Created : 20/04/2021
Time : 22:04:25
Diagnostics Date : 20/04/2021




NOTIFY PATIENT

CHF DETECTION
HF CLASSIFICATION
HEART DISEASE

DIAGNOSTIC
STATISTICS

UPLOAD

S.no	Date	Time	CHF	Heart Failure	Heart Disease
1	20/04/2021	22:04:25	Yes	-	-
2	20/04/2021	22:04:25	Yes	HFpEF	No



Home
Dashboard
Name
Logout

PA0445
Name : Vanessa Russell
Date Created : 20/04/2021
Time : 22:04:25
Diagnostics Date : 20/04/2021
NOTIFY PATIENT

CHF DETECTION

HF CLASSIFICATION

HEART DISEASE

DIAGNOSTIC

STATISTICS

UPLOAD

Diagnostics Run Successful
Diagnostics Detected : 60
Diagnostics Executed : 57

CHF	Heart Failure	Heart Disease
Yes	-	-
Yes	HFpEF	No

Home
Dashboard
Name
Logout

PA0445
Name : Vanessa Russell
Date Created : 20/04/2021
Time : 22:04:25
Diagnostics Date : 20/04/2021
NOTIFY PATIENT

CHF DETECTION

HF CLASSIFICATION

HEART DISEASE

DIAGNOSTIC

STATISTICS

UPLOAD

S.no	Date	Time	CHF	Heart Failure	Heart Disease
1	20/04/2021	22:04:25	Yes	HFpEF	No
2	20/04/2021	22:04:25	Yes	-	No
3	20/04/2021	22:04:25	Yes	HFpEF	No
4	20/04/2021	22:04:25	Yes	-	No
5	20/04/2021	22:04:25	Yes	HFpEF	No
6	20/04/2021	22:04:25	Yes	HFpEF	No
7	20/04/2021	22:04:25	Yes	HFpEF	No
8	20/04/2021	22:04:25	Yes	HFpEF	No
9	20/04/2021	22:04:25	Yes	HFpEF	No

c. Statistics (FR9)

- User can view statistics of each feature against time and against other features by clicking on the statistics button and selecting the features.

Home
Dashboard
Name
Logout

PA0445
Name : Vanessa Russell
Date Created : 20/04/2021
Time : 22:04:25
Diagnostics Date : 20/04/2021
NOTIFY PATIENT

CHF DETECTION

HF CLASSIFICATION


HEART DISEASE

DIAGNOSTIC

STATISTICS


UPLOAD


S.no	Date	Time	CHF	Heart Failure	Heart Disease
1	20/04/2021	22:04:25	Yes	HFpEF	No
2	20/04/2021	22:04:25	Yes	-	No
3	20/04/2021	22:04:25	Yes	HFpEF	No
4	20/04/2021	22:04:25	Yes	-	No
5	20/04/2021	22:04:25	Yes	HFpEF	No
6	20/04/2021	22:04:25	Yes	HFpEF	No
7	20/04/2021	22:04:25	Yes	HFpEF	No
8	20/04/2021	22:04:25	Yes	HFpEF	No
9	20/04/2021	22:04:25	Yes	HFpEF	No



d. Notify Patient (FR8)



- User can Notify Patient by clicking on Notify Patient Button and entering the email of the patient and the message and clicking on send button.





6. Profile (FR10)

- User can open their profile by clicking on their name on the navbar on top of application.
- It contains only one functionality i.e. to delete one's account.
- User can delete their account by typing "Delete" in the field and submitting it. And the User's account will be deleted and the User will be prompted with Delete Successful message.



Delete Sucessful

Your account has been deleted.